



Efficient immunohistochemical differential diagnosis of undifferentiated neoplasia


Conventional histomorphology: pleomorph and/or high-grade spindle-cell



Conventional histomorphology: epithelioid, intermediate-cell to large-cell




Conventional histomorphology: small-cell



Markers for immunohistochemical diagnostics of CUP (Carcinoma of Unknown Primary)

Marker	Predicted localisation of primary tumour	Type of carcinoma	Comments
Arginase-1	liver	HCC (Hepato Cellular Carcinoma)	more specific than HepPar1 and Glypican-3
CDX-2	colorectum, ovary, bladder (diffuse immunoreactivity)	adenocarcinomas (ovary only in case of gastrointestinal mucinous differentiation)	adenocarcinoma of stomach, oesophagus, pancreas, and biliary ducts frequently show heterogeneous immunoreactivity. Adenocarcinomas of the lung are rarely positive. (mainly mucinous carcinoma)
CK7/CK20	colorectum	adenocarcinomas (CK7-/CK20+)	
DPC4/SMAD4	pancreas, biliary ducts	adenocarcinomas	only the loss of expression is relevant for diagnosis!
GATA3	breast, urinary tract	adenocarcinomas	more sensitive for mamma as GCDFP-15 and Mammaglobin in poorly differentiated tumours, but less specific; among others mesothelioma and ductal pancreas carcinoma
GCDFP-15	breast	adenocarcinomas	high specificity but low sensitivity in poorly differentiated tumours; rare immunoreactivity in lung
Glypican-3	liver	HCC (Hepato Cellular Carcinoma)	positive in malignant melanomas, in a minority of squamous cell carcinoma, and in yolk sack tumours (> 90%); otherwise more sensitive and specific than HepPar1
HepPar1	liver	HCC (Hepato Cellular Carcinoma)	positive in appr. 50% of all adenocarcinomas of stomach and to a lesser degree in other primary localisations (colon, lung, pancreas); therefore minor positive predictive value for differentiation from HCC and liver metastases.
Mammaglobin	breast, endometrium	ductal/lobular carcinomas, endometrial adenocarcinomas	high specificity, but low sensitivity in poorly differentiated breast carcinoma
NKX3.1	prostate	adenocarcinomas	high specificity and sensitivity; positive in some PSA-neg., poorly differentiated adenocarcinomas.
Estrogen Receptor	mamma, ovary, corpus uteri	adenocarcinomas	negative in mucinous carcinomas of the ovary and serouse corpus carcinoma
Pax-8	kidney, ovary, uterus, thyroid gland	carcinoma of the kidney (all types), mullerian tumours; thyroid carcinoma (including medullary)	lower sensitivity in mucinous ovary carcinomas; more sensitive than Thyroglobulin in insular and anaplastic carcinomas of the thyroid; mesothelioma negative
SATB2	colon	adenocarcinomas	more specific than CDX-2, larger proportion of stained medullary carcinoma
PSA	prostate	adenocarcinomas	
Thyroglobulin	thyroid gland	differentiated and insular thyroid gland carcinomas	lower sensitivity in insular carcinomas
TTF-1 nuclear	lung, thyroid gland	adenocarcinomas (non mucinous), large-cell non-neuroendocrine carcinomas; insular and medullary thyroid carcinomas	with clone 8G7G3/1 also rare immunoreactivity in endometrial adenocarcinomas. expression in small cell carcinomas is not location specific
Uroplakin II	urinary tract	urothelial carcinoma	very specific, more sensitive than Uroplakin III
WT1 nuclear	ovary	serous adenocarcinoma	mesotheliomas and mucinous mamma carcinomas are positive; serous corpus carcinoma are negative

Lokalisation of the primary of neuroendocrine tumours (NET), G1 and G2			
CDX-2	midgut (ileum, appendix, colon)	NET G1/2	
NKX6.1	pancreas, duodenum	NET G1/2	alternative to polyclonal Pax-8
Pax-8 polyclonal	pancreas, duodenum	NET G1/G2	due to cross reactivity to Pax-6
TTF-1 nuclear	lung	carcinoid	



Comments: This chart shows possible algorithms for immunohistochemical differential diagnosis after conventional histomorphological analysis of largely undifferentiated adult malignant neoplasms. Certain neoplasms frequently occurring in children and adolescents (e.g., rhabdomyosarcoma) are partially not covered by the proposed marker constellations. The antibody panels are based on the assessment of tumour localization and morphology of tumor cells and allow in the majority of cases at least a basic and cost-effective classification of the lesion. In a second approach additional classificatory, prognostic and predictive markers specific to the identified histogenesis of the tumor can be used. Since none of the markers has one hundred percent sensitivity and specificity, their predictive values also depend on the relative a priori (pre-test) probabilities of the tumor entities. For example, if a tumor is rare in a certain localization, age group or gender (low a priori probability), even a relatively specific marker may have a low positive predictive value. In contrast, if the a priori probability is high, even a relatively nonspecific marker has a high predictive value for the tumor entity in question. We would like to thank Dr. med. habil. Olaf Kaufmann for his expert advice in the design of this chart. © ZYTOMED SYSTEMS 2015